



Equipment Catalogue

Contact Information:

HealthBoss/MediSCEN

Call 905-262-0927 or 905-938-0927 or toll-free 1-888-938-0927

Email support@healthboss.org

SCENAR and COSMODIC Devices

We are offering you 4 SCENAR models, which we consider "your best bet" this holiday season.

Professional SCENAR devices:

EX735Ag and EX735Ag M (Modific)

These are both The Finest SCENARs That Have It All – SCENAR-COSMODIC Hybrids containing both SCENAR and COSMODIC technologies and a software that automatically alternates them in interaction with the body.

The differences between two models are rather quantitative than qualitative. The Modific is faster and slightly more powerful, but other than that both models have similar structure and functions.

If you want the best of the best – go for the EX735Ag Modific.

If you want "the next best", which is a bit more affordable (\$1,000 difference) – go for the EX735Ag.

**And no matter what you choose, if you are a health practitioner,
the EX735 is made for you.**

Home SCENAR devices:

The D.O.V.E. Scenar and the 705m

The D.O.V.E. Scenar is our overall favourite home SCENAR device and is simply The Best Home SCENAR Ever Made. It features the purest SCENAR technology – and yet it has some elements of COSMODIC. So in a way it is a "crossbreed" between two Scenar generations and, therefore, can be used as either SCENAR or COSMODIC. You cannot go wrong with The D.O.V.E. Scenar.

The 705m is "the C-DOVE's little sister". They both feature pure COSMODIC technology. The difference between the C-DOVE and the 705m are like the difference between EX735Ag M and the EX735Ag – slight yet substantial. Here we have the same "effect vs price" relation. The C-DOVE may be more efficient, but the 705m is much more affordable as a home Scenar (\$900 difference).

As per the comparison between The D.O.V.E. Scenar and the 705m, in our opinion, The D.O.V.E. is more versatile and better positioned as a home Scenar device. The 705m is solely COSMODIC and needs to be used as such, while The D.O.V.E. Scenar is a "crossbreed" and can offer you the user benefits of both technologies.

For example, you can Dose the spot with The D.O.V.E. like you would Dose with COSMODIC – but you can also attach remote electrodes to it and do "Little Wings", the "SCENAR massage" or any other typically classical SCENAR procedure, which you cannot do with COSMODIC.

Overall, The D.O.V.E. is a somewhat "challenging" and intense Scenar device great for quick pain and stress relief, traumas and injuries, young athletic people, for dental problems, bones and joints etc.

COSMODIC devices are more delicate yet profound, and if you are dealing with internal organs disorders, auto-immune or hormonal disorders and the like – a COSMODIC unit like the 705m should work really well for you.

These differences are relative though. Both models are universal, as any Scenar device, and both will be awesome Personal Health Assistants to every family.

It's like cars: they will all take you where you need to go, but some of them will make your travel more quick and pleasurable than others.

SOLARIS Blanket

Our SOLARIS Blanket is everybody's favourite due to its universal nature and gentle yet profound healing effect on the body.

It can be used as a stand-alone therapy – and yet it enhances any therapy it is combined with.

The SOLARIS comes in 4 sizes (the Standard size being the most popular). Ask us for details.

Soft LASER – Q-Laser System

Q1000 and Q1000NG

Q1000NG (Next Generation) is the latest version of the Q1000 laser, which is now called the Legacy.

Both models have similar basic functions and the diode structure, but the NG features several upgrades (such as different power supply, Q-Cards for programming, aluminum housing etc.) Pricing is almost identical.

We are now distributing Q1000NG, but we still can get you the good old Q1000 (by special request); we have just a few left in stock.

Q10

Great "light" version of the Q1000 well suited for home use and much more affordable in price (\$2,500 difference).

660 & 808 Flash Probes

New and substantially improved versions of the old 660 & 808 Enhancers. They are still the stimulating lasers with 660 and 808 wavelengths, but new functionality and nice ergonomic design makes them look and feel almost like different devices.

They are just slightly more expensive than the old 660 & 808 Enhancers, but it is well worth it.

Web Links for More Detailed Information About the Products

- SCENARs: <http://21centurymedicinechest.com/scenar/>
- The D.O.V.E. Scenar: <http://dovescenar.com>
<http://scenar.biz/scenar-products/scenar-devices/dove-scenar.html>
- The C-DOVE/the 705m: <http://scenar.biz/scenar-products/scenar-devices/scenar-705.html>
- The EX735Ag & EX735Ag M: <http://ex735.com>
<http://scenar.biz/scenar-products/scenar-devices/scenar-cosmodic-ex735ag.html>
- Solaris Blanket: http://21centurymedicinechest.com/solaris_blanket/
- Soft Laser: <http://21centurymedicinechest.com/softlaser/>
- Q1000: <http://www.laserscenarfusion.com/products/q1000-soft-laser.php>
- Q10: <http://www.laserscenarfusion.com/products/q10-soft-laser.php>
- 660 & 808: <http://www.laserscenarfusion.com/q-series-enhancers.php>
- Q1000NG and Flash Probes: <http://softlaser.net/>

OUR SATISFACTION GUARANTEE

Your Satisfaction is 100% GUARANTEED.

Use the **Equipment** for **60 days**, if you are not convinced it's worth many times your investment you may return the equipment in its original condition and receive a full refund, less shipping and handling and a 15% restocking fee.

If you purchase a **Fusion Package** you may use the equipment for **90 days**, and if you are not convinced it's worth many times your investment you may return the equipment in its original condition and receive a full refund, less shipping and handling and a 15% restocking fee.

Like the music business individual **DVD Trainings** are not returnable, we will happily replace any damaged disk. If for any reason you are not happy with the content of a training DVD please give our offices a call, and we'll do our very best to make you happy.

Why do we charge a restocking fee? First, any of our electronic equipment returned CANNOT be re-sold after it is returned, without reconditioning. Second, we sincerely hope you use the equipment to treat yourself and your family for the full 60 - 90 days. The 15% restocking fee is significantly less than what we would rent you the equipment for over a 60-90 day period, and having it in place helps us keep the initial purchase prices lower.

HealthBoss/MediSCEN
Call 905-262-0927 or 905-938-0927 or
toll-free 1-888-938-0927
Email support@thehealingfusion.com



SCENAR and COSMODIC Devices

The D.O.V.E.

(Device Organizing Vital Energy)
Is the ultimate in PAIN MANAGEMENT
AND PERSONAL HEALTH CARE



This cell phone sized electronic wonder will **eliminate your pain in a matter of minutes** by jumpstarting and accelerating the body's natural healing processes; very much like the famous Star Trek healing machine.

The **DOVE** is the latest in the long and proven evolution of **SCENARs** - universal healing devices, originally developed within the Russian space research program in the 80's. In its essence, it is an intelligent adaptive electrical stimulator.

The **DOVE uses biofeedback** in order to establish a dialog with the body and convince it to unlock those amazing healing and regenerative powers that we all possess yet often do not use.

It is completely non-invasive; it simply enhances the body's own healing and allow it to perform those "healing miracles" you always hear of.

**The body is the best doctor around already,
and The DOVE gives it a Ph.D.**

Priority conditions for treatment:
Musculoskeletal system disorders
Digestive disorders
Dental Disorders
Disorders of the nervous system
Respiratory diseases
Emergency care

Can also be used for treatment of:
Eye disorders
Disorders of urinary system
Cardiovascular diseases
ENT disorders and much more

The DOVE is **the newest "true" SCENAR** on the market. It has been released in 2006, and quickly became one of the most popular SCENAR devices for home use.

Dr. Alexander Karasev, the inventor of SCENAR-Technology and The Master SCENAR Designer, gave the DOVE his best shot and created "**a SCENAR as it should be**" (in his own words).

After years working exclusively on COSMODIC technology, Dr. Karasev returned to the idea of designing the perfect home SCENAR and re-designed "the grandmother of all SCENARs", his original famous "troika" **SCENAR-035** giving it a much cleaner signal, some new features, a new look and a new life in The DOVE and its little Russian sister Mioscen.

The DOVE is a SCENAR with only one feedback channel – and yet it has characteristics of second-generation SCENARs in its interaction with the body. It is a sort of "crossbreed" between SCENAR and COSMODIC - and can be used equally successfully in the manner typical of both SCENAR generations. In other words, in the DOVE you have the best of both worlds, and a devoted health assistant anywhere you go: *simple, efficient, safe, and reliable.*

"Main purpose of design: Very-very simple, but universal.

Dr. Karasev on The DOVE:

It is a reincarnation of the legend. The most affordable device, which can treat many diseases. Very compact size with not too many features. Can work fully automatically or be operated manually, with the help of navigator-wheel.

People call it "a device for babushka" (there is no exact and correct English word for the vivid Russian word "babushka", so I have to use it "as is"... if you are curious what it means, welcome to Internet!). It can do a little here and a little there, when the body has numerous diseases and problems which need to be taken care of all the time. And, of course, it will help anybody, not just "babushka".

Those whose hands are not "smart" enough, can just lie on it and let it do the job.

This is a really a "people's device".

The 705m

The SMALLEST COSMODIC

Purpose:

705m is intended for **both at home and professional practice use**. It is fully automated in action and operation.



Main features:

- Switches-on on touch
- Energy level is automatic – and can be manually set.
- There is a “dynamic adaptation” mode which allows you to automatically manage stimulation during prolonged action
- You can switch the sound off and on
- Easy to determine the beneficial treatment areas

LED indication displays:

- The device contact with the skin.
- Preliminary analysis of the area – Let you know if this is a good place to treat.
- The levels of energy used
- The body’s reaction or response to the treatment
- “The Dose” when enough treatment is done
- The activation of “dynamic adaptation”
- Battery status

Priority conditions for treatment:

Emergency aid
Stress and related conditions
Children’s diseases
Eye disorders
Ear-Nose-Throat Problems
Dental Disorders

Can also be used for treatment of:

Musculoskeletal system disorders
Respiratory diseases
Digestive disorders
Disorders of the nervous system
ENT disorders
Disorders of urinary system
Cardiovascular diseases
Disorders of hematopoietic system

Dr. Karasev, the SCENAR inventor, on the 705 m:

“And now the smallest, baby SCENAR-COSMODIC. Main purpose of design: First-aid, travel.

Such a cutie! Always comes to the rescue. Very power-saving, extremely small and does not need any remote electrodes (it has two built-in electrodes which can work on both big and small areas). This makes it so easy having it always in your pocket.

*This device will always help you when you are going out, no matter what happens. Its main mission is to **relieve pain quickly** and only after that start treatment. That’s how its algorithms are arranged.*

Automatic function is extremely simple — just place it on the skin and let it do the job. Even though it is small, it is very powerful! Headache, toothache, pains in the back, bruises, cuts... or if you ate something not so good... just place it on the place that hurts and let it work. With this device you can feel secure when you work, or go on a business trip.

*You can’t help falling in love with this baby! It is unique because **it can be with you the moment you need it**. This is invaluable, when you or someone near you needs help.”*

We also have an advanced version of this model, which allows external electrodes. There is also a lithium rechargeable battery for this device, which can be charged from a charging unit supplied with it, or a cell-phone USB charger, or a USB port of your PC.

The EX735Ag & EX735Ag M

(SCENAR-COSMODIC Hybrid)

**The FINEST SCENAR
that HAS IT ALL**

**When It Comes to
a SCENAR device, the
EX735Ag is as good as it gets.**



*"I will tell you something very personal: every time I pull my EX735 out of the leather case in public and see eyes opening up in amazement at this wonder of electronic engineering, I feel proud... Proud to own **The Finest SCENAR in the World**; proud to know that in a few moments it will create yet another healing miracle; proud to be Russian and share a homeland with the device that comes as close to the Star Trek Healing Device as it is possible in the first decade of the 21st Century."*

Irina Kossovskaja, M.D., Ph.D., D.N.M.

The **EX735Ag M** combines **both SCENAR and COSMODIC** technologies and **puts therapy on autopilot** by switching between them at the right moment. Artificial intelligence in the palm of your hand.

Purpose

EX735Ag M is a top-of-the-line **professional SCENAR** device, which can also be successfully used at home. Fully automated action and operation.

New features (there is no other device which allows you achieve the following):

- automatically select healing technology** which helps you to achieve optimal effect of treatment;
- form subsequent impulses in advance, speeding up the recovery (**ACTIVE SCENAR & COSMODIC**);
- automatically regulate action energy levels in both SCENAR and COSMODIC which excludes "defence reaction" of the body;
- correct manual settings** of action modes if they are wrong;
- start treatment immediately using all functions of automatic control (turns on itself, tunes itself, etc)
- maintain **perfect fine-tuning** via user mode for automatic correction
- keep **accurate records** and regulate time of use of the device by medical personnel

Dr. Karasev, the SCENAR inventor, on The EX735Ag M:

"Main purpose of design: Treatment of variety of acute and chronic diseases on the deeper than physical levels; disorders of emotional and spiritual origin, disorders of psyche.

*This device deserves a novel! It is on top of the evolution of SCENAR history — **EX735Ag Modific!** It is the one we are the most proud of. This is the first device, in which we could unite nearly all beneficial features of other devices, and achieve results which only God could do so far. This device is our newest, but it has already gained indisputable authority.*

*Saying that it can do everything would be incorrect... SCENAR is continually evolving, and this device is not the end. However, abilities of EX735Ag Modific surprise us every day... it can even cope with Parkinson's disease and many other diseases which until now have been regarded as incurable. On the one hand, it's an **extremely sophisticated** device. On the other hand, **even a 3 year old can handle it.***

*Yes, this device is an instrument for a real professional and in skilful hands can work wonders. But when life drives you into a corner you need this device to be able to help yourself or someone else!
When every family has Modific, we all will be rich... rich in health and soul."*

The SOLARIS Energy Blanket

"I am in love with my Blanket..."

Lisa C., Mediscen Client

The **SOLARIS Blanket** has been developed almost a quarter century ago within a Russian Space Research Program as an evolution of the famous Wilhelm Reich's "organ accumulator."

It is the simplest Feedback device that represents an "energy mirror" to the body.

The SOLARIS Blanket "cocoons the body in a field of positive self-generated radiations" making you younger, stronger, happier, and, yes, **radiant**. **You have to try it to FEEL it!**

SOLARIS Energy Blanket is proven to be excellent for daily:

- Stress Relief,
- Easing of Pain and Improving Sleep,
- Prevention of Premature Aging.

SOLARIS:

- protects the body from **harmful outside energy fields**;
- reflects and **accumulates the body's own energy**;
- improves **energy distribution**;
- mirrors the body's "energetic portrait" thus providing the body with the **Feedback** it needs to regulate its own energy flow.

By gently and gradually aligning the body's energetic processes, it **assists in healing** from a wide variety of dis-ease conditions.

SOLARIS is **entirely safe and natural**. It requires no special training.

Use it for pain relief...stress relief...sleep disorders...to increase resistance to seasonal infections and viruses...to reduce the symptoms of many diseases...to increase blood circulation and relax muscles...as an anti-aging treatment...or just to feel warm and toasty.

The SOLARIS Blanket is good for **absolutely everyone** in **any circumstances**.

It is, however, especially valuable for **children** and **elderly people**, because it strengthens the body against dangers and hazards of its environment and preserves the body's energy resources. SOLARIS can be used to great advantage to maximize **athletic performance**.

It has also proven very effective for **veterinary use**.



Doc David is getting **ENERGIZED**



The Q Laser System

A Pain Free, Non-invasive Method for Cellular Healing

- ☑ Soft (Low Level) LASER Therapy works by directing a beam of **non-heating LASER light** to the area involved.
- ☑ **On a cellular level** this increases local circulation, improves tissue oxygenation and stimulates the production of ATP, allowing tissue to be more responsive and work as intended.
- ☑ On a larger scale, LASER therapy quickly decreases pain and increases the tissue's ability to handle stress.
- ☑ Low Level Laser Therapy is a completely safe, noninvasive and non-thermal therapy in which a special type of light penetrates the surface of the skin as well as underlying tissues to stimulate natural healing in the body. The body's cells are exposed to photon energy, which in turn...
 - Reduces pain
 - Reduces inflammation
 - Increases circulation
 - Increases cellular ATP production
 - Releases tight muscles (both smooth and striated)
 - Speeds bone repair
- ☑ Low Level Lasers are also used in acupuncture to stimulate the acupoints of the body without using needles.

Over the last 13 years, we have conducted a thorough research and clinical comparison of the different Soft Lasers and selected the Q-LASER System as our choice for you for the following reasons:

- **Q-LASERS (Q1000NG & Q10)** are **multi-colored**. They combine light of several wavelengths – and it is a well-known fact that different tissues in the body respond to different wavelength.
- Recipe for the color combination, which is **unique and proprietary to the Q-LASER System**, allows for the light of different wavelength not to weaken each other due to interference but to enhance and strengthen each color's effect due to forming a **SOLITON WAVE** of light in the body.
- The **SOLITON WAVE** is well-known in nature and can be seen in tidal waves, dolphin sound communication, and many other natural phenomena where a wave travels a very long distance without losing its size and strength. It is formed when several waves come together in a very particular way making a whole greater than the sum of its parts. The SOLITON WAVE of Q-LASER light travels within the body uninterruptedly carrying the nourishing energy to all organs and tissues. The Q-LASER light cascades from the place it is applied to everywhere like the body is transparent.
- The above unique quality allows the use of a **VERY LITTLE LASER POWER** making the **Q-LASERS Q1000NG and Q10 RESONATING LASERS** while most of the other Soft LASERS on the market are STIMULATING. The body likes the low LASER power and responds to it much better as it is not overwhelming to the body and nourishes without damaging. Q1000 and Q10 Resonating LASERS are so safe that you can use them successfully even directly on the eyes.
- While the Resonating LASERS are very helpful for the widest variety of disorders, some of ill-health conditions still require the power of a Stimulating LASER. **Q-LASER system** is the only one on the market today that provides for **combination of Resonating and Stimulating LASERS**. The Stimulating **660 and 808 LASER Enhances** are used in combination with Q1000 utilizing its battery power. That gives you a **unique blend of three LASERS for the price of one**.

FDA confirmed the efficiency of the Q1000NG and 660 Flash Probe Soft LASER in treatment of osteoarthritis of the hand and allows it for sale over-the-counter.

Q1000NG Soft Laser

The Q1000NG is a Soft Laser for Both Soft & Hard Tissue

It features computer-controlled energy output which stays constant and is essential to avoid the body's tendency to polarize against fluctuating energy sources.

Features:

- 20 diodes –12 true LASER diodes 8 LEDs
- 20 diodes in the Q1000NG deliver 7 different wavelengths (nanometers) simultaneously with blue, red, and infrared output
- Delivers resonating subtle energies to embrace and heal cells
- Q1000NG is the power source for our two Flash Probes
- Re-chargeable Smart Lithium Ion battery provides approximately
- 200 minutes before charging
- Size of a cordless telephone headset
- 2-year warranty standard - Lifetime warranty available

Q1000NG offers 7 different treatment modes:

- ☑ **Mode 1:** for pain, inflammation, burns, tendon, ligament, and bone injuries
- ☑ **Mode 2:** for the heart and brain
- ☑ **Mode 3:** 29 different frequencies beneficial to glands and organs
- ☑ **Mode 4:** Anti-aging and skin mode: Anti-aging, dandruff, hair treatments, scaly skin, psoriasis, scars, varicose veins, wrinkles
- ☑ **Mode 5:** Universal Healing Code: Excellent results with stroke patients, heart and vascular concerns, ADD and even was a miracle for a client with convulsions. Massage therapists will follow a mode 1 treatment with a mode 5 treatment to smooth out the edges in the energetic field
- ☑ **Mode 6:** Serotonin/Dopamine mode- Useful to gradually reduce meds. Some Q1000 users have reported help with sleeping by placing on gut or head before they go to bed
- ☑ **Mode 7:** Infection, mold and yeast: 8 frequencies known to assist. And help strengthen the immune system to deal with these issues.

Registered with the FDA as a Class I laser device.

Specifications:

Industrial strength aluminum housing milled from a solid block of aluminum. Your Q1000NG is made to last a life-time!

Safe, Class I laser device.

Patented Technology, and multiple patents pending. Patent-pending QCard system allows your Q1000NG to transform into multiple devices.

7 Healing Wavelengths. All the beneficial wavelengths Dr. Lytle uncovered during his research.

Unlimited Healing Frequencies. Patent-pending technology allows you to fully customize your Q1000NG for your needs.

Light Emitting Elements: 12 laser c

Multiple wavelengths: 470-940nm

Power per element: less than 5mW

Maximum Operational Output Power: 42mW
Treatment Area: 2.69 sq. in. (17.34 cm. sq.)

Power Density: 2.4mW/sq. cm.

Energy Output: 7.5 joules per 3-minute treatment cycle

Modulation: 20kHz max., up 63 unique frequencies
Up to 6 soliton waves



Q10 Soft Laser

Designed for home use – either by professionals as a rental device ... or by lay persons themselves, the Q10 fills a much ignored niche in the laser therapy market.

With very few exceptions, regardless of how much you pay, competing products usually only offer one wavelength of laser energy.

Why limit your chances for success?

Features:

- Four True Laser Diodes With New Cutting-Edge Soliton Wave Configurations
- Many competing products costing thousands of dollars more offer less.
- Five LED's, Including A New *Blue* Wavelength
- The Q10 is one of the only laser products in the world to combine both true laser diodes and LED's in the same delivery head. The results are superior performance.
- Five Different Wavelengths of Beneficial Energy
- With very few exceptions, regardless of how much you pay, competing products usually only offer one wavelength of laser energy. Why limit your chances for success?
- Three Unique Treatment Modes
 - Mode 1 – Pain/Inflammation
 - Mode 2 – Skin/Cosmetic Enhancement
 - Mode 3 – Universal Health Benefits
- Each mode utilizes unique frequency's to provide better results.
- Application Specific Internal Battery With *Q-Charge* Intelligent Charging System
 - Get over 300 minutes of use per battery charge. Very high quality, long-life lithium-ion battery and intelligent charger mean less problems and headaches for you.
- Bright Multi-Color Display
 - No more fumbling around trying to determine what function your laser is in.
- Sleek Ergonomic Design
 - The Q10 is designed to fit even the smallest hands comfortably.
- Extremely Portable
 - At just under 7.5" tall, 1.5" thick and only 8.2 ounces the Q10 is "take it with you anywhere" small and portable.
- Intuitive One-Button Operation
 - One button user friendly design.
- Maximum Safety – Class I, *Non-Significant Risk* Laser Device
 - No risk, no worries ... safe enough to use anywhere and on anyone
- Made In The U.S.A.



Q1000NG Flash Probes

808

The 808 Flash Probe produces 300mW of 808nm near-infrared laser light.

Studies show that this power and wavelength configuration is able to penetrate deeply and is effective in helping heal cartilage, ligaments and bone issues.

Efficacious as well for bone infections, fractures slipped and herniated discs.

Specifications:

Light Emitting Elements: Single 808nm 500mW laser diode
Operational Output Power: 300mW
Treatment Area: 0.05 sq. in. (0.3165 cm²)
Power Density: 947.87 mW/cm²
Energy Output: 54 joules per 3-minute treatment cycle
Registered as a Class IIIB laser device



660

The 660 Flash Probe send out 30mW of 660nm red laser light.

This wavelength has been shown beneficial for stimulating acupuncture points in place of needles, promoting healing in soft tissues and muscles, reducing muscle spasms, healing wounds and softening scars.

Very good for acupuncture points, scarring, in the ear, and many dental applications (in the mouth) such as gingivitis, Periodontal disease, and cold sores. Not recommended for use on the eyes.

Specifications:

Light Emitting Elements: Single 660nm 50mW laser diode
Operational Output Power: 30mW
Treatment Area: 0.05 sq. in. (0.3165 cm²)
Power Density: 110.58 mW/cm²
Energy Output: 6.6 joules per 3-minute treatment cycle
Registered as a Class IIIA laser device

